
SNAP*

*Simplicity from IDEXX**

SNAP* Beta-Lactam ST Test Kit

Validated for raw, commingled cow milk

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Not for use in NCIMS milk-regulatory programs. Not for sale in the U.S.

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IDEXX
LABORATORIES


Product and Intended Use

The SNAP* Beta-Lactam ST Test is an enzyme-linked receptor-binding assay that will detect beta-lactam antibiotics residues in raw, commingled cow milk at or below established maximum residue limits. This test will detect most beta-lactam antibiotics (see Performance Information).

Introduction

Beta-lactam antibiotics are widely used in the treatment of mastitis and other infections in dairy cattle. Failure to follow antibiotic label instructions and milk-withholding guidelines may result in antibiotic residues in milk. The SNAP Beta-Lactam ST Test is designed to simplify monitoring beta-lactam residues as part of a quality assurance program.

Kit Components

Accessories (Available through IDEXX)

- SNAPshot* DSR Reader
- SNAP* positive control

Kit Storage

Kits must be refrigerated at 2°–8°C (36°–46°F).

Sample Information

- Use raw, commingled cow milk.
- Samples must be refrigerated at 0°–10°C (32°–50°F) and tested within three days of collection.
- Thoroughly mix the sample before testing.
- DO NOT use abnormal-looking milk.

Precautions and Warnings

- Do not mix sample tubes and devices from different kits.
- Do not use kit components past their expiration dates.
- If the control spot fails to develop color, retest the sample.
- The SNAP device must be run in a horizontal position.

Before You Begin

- The SNAP* Beta-Lactam ST Test is designed for use under normal ambient conditions (15°–30°C or 59°–86°F).
- Remove the number of devices needed for the day from refrigeration and allow them to warm to room temperature for at least 15 minutes before using. Unused devices can be returned to refrigeration.
- Samples must be at 0°–10°C (32°–50°F).
- Shake the milk sample thoroughly.


Test Preparation

- Remove the SNAP device, pipette and sample tube from the bag.
NOTE: A new SNAP device will show a light blue color on the undeveloped spots in the results window and in the activation circle.
- Check that the reagent pellet is at the bottom of the sample tube. If not, tap the tube to return the pellet to the bottom.

Test Procedure

1. Pipetting the Sample

- Place the SNAP device on a level surface.
- Shake the milk sample thoroughly.
- Remove and discard the sample tube cap.
- Using the IDEXX pipette, draw up the milk sample (450 μ L \pm 50 μ L) to the indicator line.


NOTE: When pipetting, take the sample from the middle of the sample container and slowly draw the sample into the pipette to avoid air bubbles.

- Add all of the milk sample from the pipette to the tube, and then shake the sample tube to dissolve the reagent pellet.

NOTE: Do not allow the sample to remain in the tube longer than 15 seconds.

2. Adding the Sample to the SNAP Device

- Pour the contents of the sample tube into the sample well of the SNAP device and discard the tube. The sample will flow across the results window toward the blue activation circle.
- Once the sample has entered the activation circle, but before the blue circle has completely disappeared, press the activator button firmly until it is flush with the body of the SNAP device.
- After activation, allow the device to **develop for 6 minutes**, then interpret the result immediately.

NOTE: If color does not develop on the control spot, retest the sample.

Interpreting Test Results

- **SNAPshot DSR/SNAPshot Reader:** Insert the SNAP device and follow the instructions for reading the device.
- **Visual Read:** When interpreting results, position the SNAP device with the sample well to the top, as shown in the diagram. If color development does not occur in the control spot, retest the sample.


Negative Result

The sample spot is darker than the control spot.

The sample spot is equal to the control spot.

Positive Result

The sample spot is lighter than the control spot.

Performance Information

Sensitivity

The SNAP Beta-Lactam ST Test Kit has been validated to detect penicillin G at or below 4 ppb in raw, commingled cow milk.

Beta-lactam Antibiotics Detected

Penicillin, ampicillin, amoxicillin, cloxacillin, dicloxacillin, nafcillin, oxacillin, cephalonium, cephalirin, ceftiofur, cefazolin, cefoperazone and cefuroxime (not a complete list)

Sole reliance on use of this kit may result in the rejection of milk that is considered safe for human consumption. It is recommended that positive samples be further analyzed using a quantitative, drug-specific method of analysis, when available.

For technical assistance

U.S. 1-800-321-0207
Europe 00800 329 43399

*SNAP and SNAPshot are trademarks or registered trademarks of IDEXX Laboratories, Inc. and its affiliates in the United States and/or other countries.

Covered by U.S. Patent Nos.: 5,726,010; 5,726,013; 5,750,333; and/or 6,007,999. Other U.S. and/or foreign patents issued or pending.

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IDEXX Europe B.V.