Total T₄ Testing Guide: Canine Hypothyroidism

Dynamic range

 $0.5-10.0 \,\mu \text{g/dL}$ (6.4–128.7 nmol/L)

Catalyst Total T₄ results

Low	$<$ 1.0 μ g/dL	(<13.0 nmol/L)
Low Normal	$1.0-2.0~\mu\mathrm{g/dL}$	(13.0-26.0 nmol/L)
Normal	$1.0 extstyle{-}4.0~\mu extstyle{g}/ extstyle{dL}$	(13.0-51.0 nmol/L)
High	$>$ 4.0 μ g/dL	(>51.0 nmol/L)
Therapeutic	$2.1 - 5.4 \mu \mathrm{g/dL}$	(27.0-69.0 nmol/L)

Canine screening

- Dogs with low total T₄ (T₄) and evidence of nonthyroidal illness (NTI) should have the NTI addressed.
- Dogs with T₄ results in the low normal range may be hypothyroid.


In dogs with low or low normal T₄ results and with consistent clinical signs, consider one or more of the following to aid in confirming hypothyroidism:

- Free T₄ (fT₄)
- Endogenous thyroid-stimulating hormone (TSH)
- Possible thyroglobulin autoantibodies (TgAA)

Hypothyroidism therapeutic monitoring

For dogs on thyroid supplement, acceptable 4–6-hour post-pill T_4 values will generally fall in the upper end of the reference interval or slightly above.

Algorithm



Total T₄ Testing Guide: Feline Hyperthyroidism

Dynamic range

0.5–20.0 µg/dL (6.4–257.4 nmol/L)

Catalyst Total T₄ results

Subnormal $< 0.8 \ \mu \text{g/dL}$ (<10.0 nmol/L)

Normal $0.8-4.7 \mu g/dL$ (10.0-60.0 nmol/L)

Gray zone in old

or symptomatic cats $2.3-4.7 \mu g/dL$ (30.0-60.0 nmol/L)

Consistent with

hyperthyroidism $> 4.7 \mu g/dL$ (>60.0 nmol/L)

Feline screening

Cats with consistent clinical signs and total T₄ (T₄) values in the borderline high range (gray zone) may have early hyperthyroidism with a concurrent nonthyroidal illness (NTI).

In these cases, consider the following to aid in confirming the diagnosis:

- Free T₄ (fT₄)
- T₃ suppression test
- · Radionuclide thyroid imaging

Hyperthyroidism therapeutic monitoring

Following treatment with methimazole (or similar), T₄ values will generally fall within the lower to midportion of the reference interval.

Algorithm

Feline hyperthyroidism suspected Common clinical signs in cats · Weight loss **Initial database** · Hyperactivity Total T₄ Polyphagia • CBC · Palpable goiter · Chemistry with Unkempt coat electrolytes · Complete urinalysis Low T₄ Normal T₄ High T₄ $0.8 - 4.7 \, \mu \text{g/dL}$ >4.7 μ g/dL < 0.8 μ g/dL (<10.0 nmol/L)(10.0-30.0 nmol/L) (>60.0 nmol/L)Euthyroid sick Normal T₄ (Gray Zone) or iatrogenic $2.3-4.7 \mu g/dL$ (30.0-60.0 nmol/L)*If strong suspicion of hyperthyroidism still exists consider retesting in 4-6 weeks or a technetium scan. CBC = Complete blood count Low or Note: 1 µg/dL is equal to 12.87 nmol/L. A result that High fT₄ falls within the gray zone of the assay should be normal fT4 considered ambiguous. Hyperthyroidism unlikely* Hyperthyroidism likely

