IDEXX Cornerstone* Getting Started Guide





Proprietary Rights Notice

Information in this document is subject to change without notice. Companies, names and data used in examples are fictitious unless otherwise noted. No part of this document may be reproduced or transmitted in any form or by any means, electronic, mechanical or otherwise, for any purpose, without the express written permission of IDEXX Laboratories. IDEXX Laboratories may have patents or pending patent applications, trademarks, copyrights or other intellectual or industrial property rights covering this document or subject matter in this document. The furnishing of this document does not give a license to these property rights except as expressly provided in any written license agreement from IDEXX Laboratories.

© 2014 IDEXX Laboratories, Inc. All rights reserved. • 06-28486-05

*Cornerstone, IDEXX SmartOrder, IDEXX SmartService, IDEXX Monitor, IDEXX VetLab, IDEXX VetVault, Image Coach, I-Vision Mobile, LaserCyte, Patient Clipboard, Pet Health Network, SmartLink, VetConnect and VetTest are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries.

DICOM is the registered trademark of the National Electrical Manufacturers Association for its standards publications relating to digital communication of medical information. Microsoft, Hotmail, Internet Explorer, Outlook, Windows, Windows Server and Windows Vista are trademarks or registered trademarks of Microsoft Corporation in the United States and/or other countries. Accelerated Payment Technologies and XCharge are trademarks or registered trademarks of Accelerated Payment Technologies, Inc. Topaz, ClipGem, SigLite and SignatureGem are trademarks of Topaz Systems, Inc. in the U.S. and/or other countries. AMD Athlon, Opteron and Sempron are trademarks of Advanced Micro Devices, Inc. Celeron, Intel, Intel Atom, Pentium and Xeon are trademarks of Intel Corporation in the U.S. and/or other countries. Dell is a trademark of Dell Inc. Abaxis and VetScan are registered trademarks of Abaxis, Inc. LifeLearn and ClientEd Online are trademarks or registered trademarks of VetLearn, Inc. Gmail is a trademark of Google, Inc. MWI is a registered trademark of MWI Veterinary Supply, Inc. Antech and Antech Diagnostics are registered trademarks of VCA Antech, Inc. Backup Exec is a trademark or registered trademark of Symantec Corporation or its affiliates in the U.S. and other countries. Backup MyPC is a trademark or registered trademark of Corel Corporation and/or its subsidiaries in Canada, the United States and/or other countries. TapeWare is a registered trademark of Barracuda Networks, Inc.

All other product and company names and logos are trademarks or registered trademarks of their respective holders.

IDEXX Computer Systems One IDEXX Drive Westbrook, Maine 04092 USA idexx.com/cornerstone

Contents

Click within the contents listing below to go directly to a section or page.

Introduction	5
Documents to Help You Get Started with Cornerstone 8.3	5
1: Hardware and Operating System Guidelines	7
General Recommendations	7
Operating System Compatibility	7
Hardware and Operating System Recommendations	8
Peer-to-Peer Servers and Stand-alone Computers	8
Dedicated Servers	8
Thick-Client Workstations (Including Laptop and Tablet Computers)	9
Thin-Client Workstations	9
Evaluating Your Current System Configuration	9
Tracking Your Hardware and Software Performance	9
2: Prerequisites for New Features	
Client Communications	10
Laboratory	12
Imaging 14	
Credit Card Processing	14
Inventory 15	
3. Requirements if Installing from Physical Media (DVD)	
Installing Across the Network from a Shared DVD Drive	16
4: System Preparation	
Installation Notes	18
Laboratory	18
Imaging 18	
Disabling the Windows [®] Firewall	19
Disabling User Account Control	
Changing to Install Mode (Windows Server 2008 R2)	
Other Steps to Take Before You Begin the Installation	
5: Installation on a Server or Stand-alone Computer	
Installing Cornerstone 8.3 from a DVD	
Upgrading to Cornerstone 8.3 from the Internet using IDEXX SmartService* Solutions	25
6: Installation on Workstations	
Upgrades from Cornerstone 7.6, 8.0, 8.1, or 8.2	
Installing on Workstations Using Install from Server	
Installing on Workstations from a Shared Network Folder	29
New Workstations or Upgrades from Cornerstone 7.5 and Earlier Versions	29
7: Installation for the Mobile Computing Module	
Installing Cornerstone 8.3 for the Mobile Computing Module	

If You Cannot Bring Each Laptop/Tablet into the Practice	
8: Security and Setup	35
Suggested Security by Staff and Role for New Windows in Cornerstone 8.3	35
9: Communications and Training	
Training Checklist	
Training Files	
Appendix: Starting the Database	41
Starting the Database	41
Automatic Start	41
Manual Start	41

Introduction

Welcome to the new release of the IDEXX Cornerstone* 8.3 Practice Management System. This release provides a number of user-requested enhancements, including features that provide faster access to information and a "smarter" patient treatment ordering process.

Begin implementing Cornerstone 8.3 by following this plan, which provides high-level steps for

preinstallation, installation, and postinstallation. The



benefit of following these steps is that you'll be prepared to set up and use the features, prepare your staff, avoid technical issues, and utilize the power of this release.

Documents to Help You Get Started with Cornerstone 8.3

To help ensure a smooth transition to your new IDEXX Cornerstone 8.3 software, we recommend using the *Cornerstone 8.3 Getting Started Guide* in combination with the *What's New in Cornerstone 8.3* document. Use the sections of the two documents in the order listed below.

Plan		
1	Getting Started Guide: Preinstallation (sections 1–3)	 Review hardware/operating system guidelines. Review a checklist of preinstallation tasks to prepare for and take advantage of the new features in this release. Review special requirements if you are using a DVD to install Cornerstone.
Install	and set up	·
2	Getting Started Guide: Installation (sections 4–7)	Follow the recommended steps for installation.
3	What's New in Cornerstone 8.3: Setup instructions at the beginning of each feature section	View detailed information about setting up and using the new features and functionality in this release.
4	Getting Started Guide: Security and Setup (section 8)	Reference this plan to set up security for new features and to configure of default settings.
Use n	ew features	
5	Getting Started Guide: Communications and Training (section 9)	Introduce the new features you choose to implement to your practice, and train staff.
6	What's New in Cornerstone 8.3: Use instructions within each feature section	View detailed information about how to use the new features.

To access these resources, visit the Cornerstone Resource Library online at

idexx.com/cornerstoneresources. Copies of both documents are also installed on your computer along with Cornerstone.

If upgrading from previous versions other than Cornerstone 8.2 and to receive guidance on an upgrade plan, please call your Cornerstone Sales Representative at 1-800-283-8386 for information.

1: Hardware and Operating System Guidelines

For optimal performance, please ensure that your system meets the recommended specifications provided in this chapter. This will ensure the best possible operating speed and level of efficiency when using the IDEXX Cornerstone* Practice Management System.

Note: Cornerstone software may function when using a computer with less than the recommended specifications. If concerned that your current system may not operate at expected levels of performance, please contact IDEXX Cornerstone Support at 1-800-695-2877.

Business Class Computers, Printers and Peripherals	IDEXX recommends business class computers, printers and peripherals to ensure the best performance, reliability and warranty. Business class devices are available from manufacturers like Dell [®] and include components, configuration options and a warranty designed to meet the needs of a demanding business environment.
Network	IDEXX recommends computers and switches capable of gigabit speed with a wired connection. (A wired connection will always provide a faster and more stable Cornerstone experience.) A wireless connection should be used only for mobile computers.
Antivirus (Malware Protection)	IDEXX recommends antivirus software that does not use a software firewall (or offers the option of turning off the software firewall). For security purposes, IDEXX recommends a hardware firewall (included with most routers) instead of a software firewall.
Monitor Resolution	IDEXX recommends 1024×768 or higher resolution for optimal viewing of the Cornerstone software.
Cornerstone Backup Verification Assistant	The Cornerstone Backup Verification Assistant requires one of the following types of backup software: Backup4all, Backup Exec [™] , BackUp MyPC [™] , PC Backup, TapeWare [®] or IDEXX VetVault* Backup Solution.

General Recommendations

Operating System Compatibility

Windows 98, Windows NT, Windows 2000 and Windows 2000 Server	Cornerstone software will NOT install on a computer that uses a Windows [®] 98, Windows [®] NT, Windows [®] 2000, or Windows [®] 2000 Server operating system. If your practice still uses one of these operating systems, please call your IDEXX Practice Solutions Representative for Supplies at 1-888-224-4408 to discuss your upgrade options.
Windows [®] XP 32-bit and Windows Server [®] 2003 32-bit	 In Cornerstone 8.1 or later, the ability to create patient CDs, archive images to a CD, or export cases to CD (Diagnostic Imaging module only) requires Microsoft[®] .NET Framework 3.5 Service Pack 1, which is available in Windows[®] XP Service Pack 2 or Windows Server[®] 2003 Service Pack 2 or later operating systems. Windows XP and Windows Server 2003 can be upgraded to the latest Windows service pack through Windows Update. With the latest Windows service pack installed, you can then download Microsoft .NET Framework 3.5 Service Pack 1.
Windows Vista Business 32-bit and Ultimate 32-bit	Windows Vista [®] requires Cornerstone 7.6 or later for workstations and Cornerstone 8.0 or later for peer-to-peer servers.

Windows 7 Professional 32/64-bit and Windows Server 2008 R2	Windows [®] 7 Professional and Windows Server [®] 2008 R2 require Cornerstone 8.0 Service Pack 2 or later. IDEXX supports both 32-bit and 64-bit versions of Windows 7 Professional but recommends 64-bit because of memory restrictions and driver compatibility with Windows Server 2008 R2, which is 64-bit only. (Note: Some compatibility issues have been identified for using certain third-party applications and devices with these operating systems. A PDF document detailing these issues is available at <u>idexx.com/cornerstonespecifications</u> .)
Windows Internet Explorer	Windows [®] Internet Explorer [®] 7 or later must be installed on each workstation. It does not have to be the default Internet browser.

Hardware and Operating System Recommendations

Peer-to-Peer Servers and Stand-alone Computers

Peer-to-peer servers allow up to 6 thick-client workstations and run a non-server operating system.

Operating System	° Windows [®] XP Professional 32-bit
	 Windows[®] 7 Professional 64-bit
Processor	Multicore Intel [®] processor or AMD processor [†]
Hard Drive	200 GB
Memory	2 GB
	° DVD-RW
Backup Device	 Removable disk drive (Dell RD1000)
	 IDEXX VetVault* Backup Solution

[†]Excluding Intel[®] Celeron[®], Intel Atom[™] and AMD Sempron[®] processors.

Dedicated Servers

Dedicated servers allow for 7 or more workstations and run a server operating system. This type of server should not be used as a Cornerstone workstation.

Operating System	 Windows Server 2003 32-bit 	
	 Windows Server 2008 R2 Standard/Enterprise 	
Processor	Multicore Intel [®] Xeon [®] or AMD Opteron [®]	
Hard Drive	250 GB [‡]	
Memory	4 GB ^{\$}	
	◦ Tape backup (LTO)	
Backup Device	 Removable disk drive (Dell RD1000) 	
	 IDEXX VetVault Backup Solution 	

[‡]Servers that will be storing digital images may require more hard drive and backup space depending on usage.

[§]Servers that will be running Remote Desktop Services (formerly known as Terminal Services) will require more memory. An additional 400–500 MB per session is recommended, depending on usage and additional software.

Thick-Client Workstations (Including Laptop and Tablet Computers)

Operating System	 Windows XP Professional 32-bit Windows 7 Professional 64-bit
Processor	Multicore Intel processor or AMD processor ^{††}
Hard Drive	60 GB
Memory	2 GB

Thick-clients are workstations that run Cornerstone locally.

^{††}Excluding Intel Celeron, Intel Atom and AMD Sempron processors.

Thin-Client Workstations

Thin-clients, also known as terminals, are workstations that run Cornerstone on the server using a Remote Desktop Connection. Thin-clients are not recommended for viewing digital images.

	 Windows XP Embedded
Operating System	 Windows XP Home
	 Windows 7 Home Premium
Processor	Intel Celeron, Intel Atom or AMD Sempron
Memory	512 MB

Evaluating Your Current System Configuration

If you are using a Dell computer system, we provide access to Dell's online tool to quickly analyze your system and help you determine if your system meets the hardware requirements for this release of Cornerstone software.

To analyze your system:

- 1. Access the Dell[®] System Analyzer at <u>accessories.us.dell.com/sna/DellPartsFamily.aspx</u>.
- 2. Click the Analyze System button to start the process.

Tracking Your Hardware and Software Performance

The IDEXX Monitor* Cornerstone* Dashboard report tracks hardware and software performance to help practices ensure that Cornerstone operates smoothly and efficiently. This monthly report is available at no charge to IDEXX Cornerstone users who have enabled the data services agreement in Cornerstone (on the **Tools** menu, select **Data Services Agreement**).

To register for the IDEXX Monitor Cornerstone Dashboard report or for more information, call 1-888-557-6517.

Many of the new features in the 8.3 release require some setup in the Cornerstone* Practice Management System before you can use them. In addition, a few features have more involved prerequisites, such as purchasing an additional Cornerstone module, or activating an account to use other IDEXX or non-IDEXX software or devices. The table below gives an overview of new features with significant prerequisites or setup requirements.

If you have questions or would like more information, please contact your IDEXX Cornerstone Sales Representative at 1-800-283-8386.

For more detailed information about how to set up and use the new features, see *What's New in Cornerstone 8.3.* Places to find this information are listed in the "For full details" sections of the tables below.

Note: Modules marked **Reinstall** are activated by reinstalling the Cornerstone software with a new activation key. To save time, you may want to purchase these modules before upgrading to Cornerstone 8.3. You can enter the new key as part of the upgrade process, rather than reinstalling again later.

Client Communications

IDEXX Pet Health Network* Pro

IDEXX Pet Health Network* Pro lets you manage communication to and from your clients within the Cornerstone* Practice Management System. IDEXX Pet Health Network Pro is currently available in the United States only.

The tools available when you use IDEXX Pet Health Network Pro with Cornerstone include:

- Sending automatic reminders by mail and/or email
- Allowing clients to confirm via email or text-with the confirmation recorded in Cornerstone automatically
- Manage important client communications from the Cornerstone message center—including appointment
 and prescription requests and updated client information
- Sharing selected documents from the patient record with clients electronically

Prerequisites

Setup

Before you use the IDEXX Pet Health Network Pro □ When you sign up for IDEXX Pet Health Network Pro, your dedicated communications specialist features in Cornerstone, you need to sign up for the service. Get started by contacting an works with your practice to customize your IDEXX Pet Health Network Pro communications marketing plan and ensure that you get the most out specialist at 1-855-627-6570 or of your client communications. The specialist will support@pethealthnetwork.com. also walk you through setting up IDEXX Pet Health Network Pro to support your plan for your practice. Note: IDEXX Pet Health Network Pro is a subscription-based service. There is also a □ If you are already using IDEXX Pet Health Network charge for postcard communications. Contact Pro without Cornerstone, some additional setup is IDEXX Pet Health Network Pro for subscription required to use the features that are available with information. Cornerstone. Contact an IDEXX Pet Health Network Pro communications specialist at 1-855-627-6570 or □ To communicate with IDEXX Pet Health Network support@pethealthnetwork.com. Pro, the Cornerstone server must have Internet access.

For full details

See "IDEXX Pet Health Network Pro" in What's New in Cornerstone 8.3.

IDEXX Email Service and IDEXX Faxing Service

Cornerstone has enhanced its email capability and added the ability to fax documents directly from within the system, saving staff members' time and enhancing the overall work flow of the practice. This feature may be of particular interest to specialty/referral and emergency practices, as well as to referring veterinarians.

Prerequisites	Setup
 To purchase the IDEXX Faxing Service, contact IDEXX Cornerstone Support at 1-800-695-2877. The IDEXX Faxing Service is available in the United States and Canada. 	To activate IDEXX Faxing Service, you will need to perform a one-time setup to specify your account information and to agree to the terms and conditions of use.
 No special prerequisites are required to set up IDEXX Email Service. As with any Cornerstone email feature, you need 	To activate IDEXX Email Service, you may need to perform a one-time setup to specify your account information.
to have an email account. If you want Cornerstone to use a desktop email program, such as Microsoft [®] Outlook [®] , that program needs to be set up.	If you want to use a web-based email provider, such as Gmail [™] , Hotmail [®] , etc., additional configuration may be needed.

For full details

See "Communications Management for Email and Fax," in What's New in Cornerstone 8.3.

LifeLearn[®] ClientEd Online™

Cornerstone integrates with LifeLearn[®] ClientEd Online[™] to offer customizable client education documents. ClientEd Online has more than 1,200 documents in its library, and new documents are added regularly. LifeLearn ClientEd Online is available in the United States and Canada.

Prerequisites

- Both new and existing LifeLearn customers will need to subscribe to the ClientEd Online service.
 Please contact your IDEXX Cornerstone Sales
 Representative at 1-800-283-8386 for more information.
- □ Using LifeLearn ClientEd Online documents requires Microsoft® Word 2003 or newer.

Setup

- Once you have completed the signup process and received your LifeLearn username and password, you can then perform the synchronization in Cornerstone to begin printing your ClientEd Online client handout documents.
- If you are currently using CD-based LifeLearn documents, you will want to review and potentially delete these documents before performing the synchronization. See "Getting Started with LifeLearn ClientEd Online" in What's New in Cornerstone 8.3.

For full details

For information about purchasing LifeLearn ClientEd Online or about integrating your ClientEd Online documents with Cornerstone, see "Getting Started with LifeLearn ClientEd Online" in *What's New in Cornerstone* 8.3.

For full details, see "LifeLearn ClientEd Online Integration with Cornerstone" in What's New in Cornerstone 8.3.

Laboratory

VetConnect* PLUS

VetConnect* PLUS is an exciting, new web-based tool that lets you view, trend, and compare current and historical IDEXX test results side-by-side. With VetConnect PLUS, you can trend in-house diagnostics and outside laboratory data together, and display the results relative to the reference intervals for both, providing an enhanced diagnostic view for you and for your clients.

VetConnect PLUS is integrated into the laboratory work flow you use in Cornerstone—giving you easy access when and where you need it. VetConnect PLUS is currently available in the United States only.

Prerequisites

Setup

- To see results from an IDEXX VetLab* Station on VetConnect PLUS, IDEXX SmartLink* technology is required to provide bidirectional flow between Cornerstone and your IDEXX VetLab Station. In addition, IDEXX SmartService* Solutions needs to be set up for the IDEXX VetLab Station. For more information, please contact your IDEXX Connectivity Consultant at 1-877-256-8457.
- □ The workstation where you want to view VetConnect PLUS must have Internet access.
- Windows[®]Internet Explorer[®] 7 or later must be installed on each workstation. It does not have to be the default Internet browser.
- You must have a VetConnect account with administrator user name and password. Visit vetconnect.com or call 888-433-9987 to request an account if your practice does not already have one.

- Activate VetConnect PLUS from within Cornerstone. This begins the process of loading historical IDEXX Reference Laboratories data (if applicable); this may continue for up to 24 hours. (Cornerstone performance is not affected by the upload.)
- After you activate VetConnect PLUS from within Cornerstone, activate VetConnect PLUS from your IDEXX VetLab Station (if applicable). This begins the process of loading data, which may take a few hours. You can still use the IDEXX VetLab Station and Cornerstone during the upload.

For full details

See "VetConnect PLUS" in What's New in Cornerstone 8.3.

IDEXX VetLab* Station

Cornerstone 8.3 offers more convenient communication between the workstations of your Cornerstone network and one or more IDEXX VetLab Station computers.

Prerequisite:

Required Setup for All Cornerstone/IDEXX VetLab Station Users

□ If you use IDEXX SmartLink technology to connect Cornerstone and IDEXX VetLab Station, there is a required one-time update to link some patient sex and species values in the database. Cornerstone will attempt to complete the update automatically, but you may be prompted to manually resolve some linking questions. See "Linking Cornerstone Values with Laboratory Values" in *What's New in Cornerstone 8.3*.

Setup for Enhanced Communication Features

 □ If you have multiple IDEXX VetLab Station computers, you can associate specific practice locations or workstations to use with a specific IDEXX VetLab Station by default. See "Improved Communications with IDEXX VetLab Station" in What's New in Cornerstone 8.3.

For full details

See "IDEXX VetLab Station" in What's New in Cornerstone 8.3.

Antech Diagnostics® Module

The Antech Diagnostics[®] module has been updated to enable online test ordering and automatic result download. In addition, pricing updates are pushed through Cornerstone automatically, providing you with the most up-to-date prices for tests. The Antech Diagnostics module is available in the United States and Canada.

Prerequisites

- To purchase the Antech module or to upgrade from the previous version of the module, contact your IDEXX Cornerstone Sales Representative at 1-800-283-8386.
- Reinstall: If you do not already have the Reference Laboratory module, you will need to purchase it as well. You will need to reinstall Cornerstone in order to start using the Reference Laboratory module.

Setup

- □ A one-time setup process is required to enter login information for the Antech website.
- There is a required one-time update to link some patient sex and species values between Cornerstone and Antech. See "Linking Cornerstone Values with Laboratory Values" in What's New in Cornerstone 8.3.
- One-time setup is required if you want to set or change default pricing options on the Antech Diagnostics Lab Information window, or if you want to manually set profile prices or change the markup percentage.

For full details

See "Antech Diagnostics Module" in What's New in Cornerstone 8.3.

Imaging

DICOM[®] Modality Worklist Integration

Cornerstone 8.3 leverages the modality worklist features of the DICOM[®] standard to share patient and exam information with a variety of non-IDEXX imaging devices. With devices that support the modality worklist and modality performed procedure step (MPPS) features of DICOM, you can send requests for images to these devices and see the status of these requests in Cornerstone. Limited integration may be possible with devices that do not support MPPS.

Prerequisites	Setup
To use the new imaging integration features, DICOM Imaging and the Modality Worklist Integration module are required. For more information, contact your IDEXX Cornerstone Sales Representative at 1-800-283-8386.	□ A one-time setup is required to enter DICOM information about your computer network and each validated system with which you want to establish a modality worklist connection. There is a one-time fee for setting up the integration between a validated
Reinstall: If you purchase DICOM Imaging, you will need to reinstall Cornerstone.	system and Cornerstone. For more information, contact your IDEXX Cornerstone Sales Representative at 1-800-283-8386. For a list of
If you have an IDEXX Digital Imaging system and already have IDEXX SmartLink Digital Imaging, you do not need to purchase additional modules.	systems that have been validated by IDEXX as capable of integrating with Cornerstone, see idexx.com/validatedsystems.

For full details

See "Worklist Integration" in What's New in Cornerstone 8.3.

Credit Card Processing

Integrated Payment Solution with XCharge®

The IDEXX Cornerstone* integrated payment solution, introduced in Cornerstone 8.1, provides a simple and secure solution for credit card processing initiated from Cornerstone. XCharge® software, developed by Accelerated Payment Technologies®, processes both credit and debit card transactions in the United States and Canada through an Internet connection. Helping you move toward a chartless practice, Cornerstone 8.3 now allows you to capture a client's electronic signature for all credit card payments using a supported signature capture device or tablet PC.

Prerequisites	Setup
 Supported signature capture devices include the Topaz[®] ClipGem[®], Topaz SignatureGem[®] (SigGem) and Topaz SigLite[®] LCD devices, or a tablet PC. For information on purchasing a signature capture device, call your IDEXX Practice Solutions Representative for Supplies at 1-888-224-4408. 	A one-time setup process is required to connect and enable the signature capture device.

For full details

See "Electronic Signature Capture with XCharge" in What's New in Cornerstone 8.3.

Inventory

IDEXX SmartOrder* Inventory Solution

IDEXX SmartOrder* is the automated inventory integration between Cornerstone and MWI® Veterinary Supply (mwivet.com). IDEXX SmartOrder allows practices to place orders within Cornerstone and receive real-time inventory item status updates and electronic record of receipt from mwivet.com via a cloud-based Internet connection. This IDEXX SmartOrder inventory solution is available in the United States.

Prerequisites

Setup

- □ An account with MWI is required to set up and use the IDEXX SmartOrder solution. Contact your MWI your practice does not have one.
 - New MWI customers must contact MWI to set • up your account and activate online ordering with mwivet.com. If you do not have an account set up, contact your MWI representative.
 - Existing MWI customers can use their current MWI customer number to activate electronic inventory ordering in Cornerstone, but you must also contact MWI to enable online ordering through MWI.
- □ A consistent, high-speed Internet connection is required for IDEXX SmartOrder.
- □ Windows[®]Internet Explorer[®] 7 or later must be installed on each workstation. It does not have to be the default Internet browser.

representative to set up your mwivet.com account if

□ To connect IDEXX SmartOrder with the MWI online inventory ordering system, you must complete a one-time setup to enter your account information in Cornerstone and to agree to the

terms and conditions of use for MWI.

□ For items that you order through MWI, you need to add a vendor item ID numbers and assign them to corresponding inventory item ID numbers. You can do this "on the fly" while working with an inventory order, or you can add and assign vendor item ID numbers ahead of time. For convenience, you may want to assign IDs for frequently-ordered items when you initially set up IDEXX SmartOrder.

For full details

See "Electronic Inventory Ordering with IDEXX SmartOrder and MWI Veterinary Supply" in What's New in Cornerstone 8.3.

3. Requirements if Installing from Physical Media (DVD)

The installation process is different if you use a DVD to install or upgrade the Cornerstone* Practice Management System.

If you have downloaded the Cornerstone* Practice Management System over the Internet using IDEXX SmartService* Solutions, skip this section and go to chapter 5, "Installation on a Server or Stand-alone Computer," starting on page 21.

If you are installing from DVD and all the computers on which you want to install the new version of Cornerstone have DVD drives, skip this section and go to chapter 5, "Installation on a Server or Standalone Computer," starting on page 21.

If you are installing from DVD, but not all of your computers have DVD drives, and:

- If a DVD drive is available on the server, be sure to install Cornerstone 8.3 on the server first, following the instructions in chapter 5, "Installation on a Server or Stand-alone Computer."
- If a DVD drive is not available on the server, you can install across the network from a shared DVD drive on another computer. See the "Installing Across the Network from a Shared DVD Drive" section that follows for detailed instructions.

If you do not have a computer with a DVD drive, please contact your IDEXX Practice Solutions Representative for Supplies at 1-888-224-4408 upgrade or alternate media solutions.

Cornerstone software will NOT install on a computer that uses any of the following older operating systems: Windows[®] 98, Windows[®] NT, Windows[®] 2000, or Windows[®] 2000 Server. For information on compatible operating systems, see section 1, "Hardware and Operating System Guidelines," on page 7.

Installing Across the Network from a Shared DVD Drive

If a DVD drive is not available on the server, you can install across the network from a shared DVD drive on another computer.

To share a DVD drive (Windows[®] XP or Windows[®] Server 2003):

- 1. On the Windows[®] XP or Windows Server 2003 computer with the DVD drive, right-click **Start** and select **Explore**.
- 2. In the navigation pane on the left, locate and click the **My Computer** icon.
- 3. In the pane on the right, right-click the DVD drive and select **Properties**.
- 4. Click the **Sharing** tab.
- On the Sharing tab, select the Share this folder option and note the Share name for this folder (typically D) for later reference.
- 6. Click OK.
- 7. In the navigation pane on the left, right-click the My Computer icon and select Properties.
- 8. Click the **Computer Name** tab and note the full computer name for later reference. (Note: The computer name is listed before the first period.)
- 9. Click Cancel.

To share a DVD drive (Windows[®] 7 or Windows[®] Server 2008 R2):

- 1. On the Windows[®] 7 or Windows[®] Server 2008 R2 computer with the DVD drive, right-click **Start** and select **Open Windows Explorer**.
- 2. In the navigation pane on the left, locate and click the **Computer** icon.
- 3. In the pane on the right, right-click the DVD drive and select Properties.
- 4. Click the **Sharing** tab.
- 5. On the **Sharing** tab, click the **Advanced Sharing** button.
- 6. Click the **Share this folder** check box and note the share name for this folder (typically D) for later reference.
- 7. Click OK.
- 8. In the navigation pane on the left, right-click the Computer icon and select Properties.
- 9. Note the computer name for later reference.
- 10. Click Cancel.

To install across the network from the shared DVD drive (all operating systems):

- 1. On the computer with the shared DVD drive, insert the Cornerstone 8.3 installation DVD.
- 2. On the server, right-click Start and then select Explore or Open Windows Explorer.
- In the Address bar near the top of the window, type \\computername\D\setup.exe (replace "computername" with the full name of the computer that has the DVD drive, and if necessary, replace D with the share name of the DVD drive). The Welcome window appears.
- 4. From this point, you can follow the standard procedure for installing Cornerstone on a server or stand-alone computer, which starts on page 21.
- 5. After the installation is complete, do the following on the computer with the DVD drive:
 - a. Remove the Cornerstone 8.3 installation DVD from the DVD drive.
 - b. If desired, turn off sharing of the DVD drive by accessing the DVD drive's **Sharing** tab and selecting the **Do not share this folder** option (for Windows XP or Windows Server 2003) or deselecting the **Share this folder** check box (Windows 7 or Windows Server 2008 R2).

4: System Preparation

Installation Notes

- Installation may take from several minutes to over an hour, depending on your computer configuration, database size and current version of the Cornerstone* Practice Management System.
- If you are upgrading from an **older version prior to Cornerstone 5.6**, call IDEXX Cornerstone Support (1-800-695-2877) to receive alternate upgrade instructions.
- If you are using the **Mobile Computing module**, see "7: Installation for the Mobile Computing Module," on page 33 in this guide.
- End of Day processing may take longer than usual the first time it is run following the installation of your new release of Cornerstone software.

Laboratory

- If your practice uses the **IDEXX SmartLink* In-house Laboratory module**, you must shut down the IDEXX VetLab* Station and all other connected analyzers before installation.
- On the IDEXX VetLab Station 3.10 or later, in the Practice Management tab's *Transmit Results* area, select the On (Transmit result records and report) option. This will set the IDEXX VetLab Station to transmit both result records and their associated IDEXX Lab Report to the Cornerstone software. (See the *What's New in Cornerstone 8.3* for more information.)
- For practices that use the following in-house laboratory devices with a serial connection, follow the applicable steps before installation: (These instructions do not apply to practices that have the IDEXX VetLab Station with IDEXX SmartLink* technology.)
 - IDEXX LaserCyte* Hematology Analyzer— Close Vetcom32 on the Cornerstone Lab workstation (CTRL+ALT+DELETE to start Task Manager > Processes tab > select Vetcom32.exe > click End Process).
 - IDEXX VetTest* Chemistry Analyzer— Close Vetcom32 on the Cornerstone Lab workstation (CTRL+ALT+DELETE to start Task Manager > Processes tab > select Vetcom32.exe > click End Process).
 - Abaxis[®] VetScan[®] Analyzer— Close Vsccom32 on the Cornerstone Lab workstation (CTRL+ALT+DELETE to start Task Manager > Processes tab > select Vsccom32.exe > click End Process).

Imaging

- If your practice has the **IDEXX Digital Radiography module**, please do the following on the workstation that is connected directly to your IDEXX Digital Imaging system:
 - 1. Open the lib folder within your Cornerstone directory (e.g., C:\Cstone\lib).
 - 2. Make a copy of the **CaptureWorkflow.properties** file. (Right-click the file, select **Copy**, then go to your Windows desktop, right-click and select **Paste**.)
- The ability to create patient CDs, archive images to a CD or export cases to CD (Diagnostic Imaging module only) requires Microsoft[®].NET Framework 3.5 Service Pack 1, which is available in Windows[®] XP Service Pack 2 or Windows Server[®] 2003 Service Pack 2 or later operating systems.

If necessary, Windows XP and Windows Server 2003 can be upgraded to the latest service pack through Windows Update. With the latest Windows service pack installed, you can then download <u>Microsoft .NET Framework 3.5 Service Pack 1</u> (also available at idexx.com/microsoftnet).

Disabling the Windows[®] Firewall

Before installation, verify that the Windows Firewall has been disabled.

Notes:

- A hardware firewall (router) is recommended in place of a software firewall.
- The following instructions assume the Windows Start menu and Control Panel are set to Classic View/Small Icons mode.

To disable the Windows Firewall in Windows XP and Windows Server 2003:

- 1. Click Start > Settings > Control Panel.
- 2. Double-click Windows Firewall.
- 3. Select Off (not recommended).
- 4. Click **OK** and close the Control Panel.
- 5. Disable any other firewall software.

To disable the Windows Firewall in Windows Vista[®]:

- 1. Click Start > Settings > Control Panel.
- 2. Double-click Windows Firewall.
- 3. Click Change Settings.
- 4. In the Windows Firewall Settings window, click Off (not recommended) and click OK.
- 5. Close the Windows Firewall window and close the Control Panel.
- 6. Disable any other firewall software.

To disable the Windows Firewall in Windows[®] 7 and Windows Server[®] 2008 R2:

- 1. Click Start > Control Panel and then click Windows Firewall.
- 2. Click Turn Windows Firewall on or off.
- 3. In the Home or Work (Private) network location settings area and/or the Domain network location settings area, click **Turn Off Windows Firewall (not recommended)**.

Note: If the options are grayed out, ensure that Turn Off Windows Firewall (not recommended) is selected and continue to the next section.

- 4. Click **OK**.
- 5. Close the Windows Firewall window and close the Control Panel.
- 6. Disable any other firewall software.

Disabling User Account Control

For Windows Vista, Windows 7, and Windows Server 2008 R2 operating systems, before installation, verify that User Account Control has been disabled.

Note: The following instructions assume the Windows Start menu and Control Panel are set to Classic View/Small Icons mode.

To disable User Account Control in Windows Vista:

- 1. Click Start > Settings > Control Panel.
- 2. Double-click User Accounts.
- 3. Click Turn User Account Control on or off.

- 4. Click **Continue** at the User Account Control prompt.
- 5. Uncheck Use User Account Control (UAC) to help protect your computer.
- 6. Click **OK** and restart the computer.

To disable User Account Control in Windows 7 and Windows Server 2008 R2:

- 1. Click Start > Control Panel and then click User Accounts.
- 2. Click Change User Account Control settings.
- 3. Slide the adjustment bar to Never notify.
- 4. Click **Yes** at the User Account Control prompt.
- 5. Click OK and click Yes! Restart the computer; otherwise, restart the computer manually.

Changing to Install Mode (Windows Server 2008 R2)

In Microsoft[®] Windows Server[®] operating systems, **install mode** is used to install applications and software, and **execute mode** is used to run those applications and software. If you are installing on Windows Server 2008 R2, you must change to install mode before installing Cornerstone 8.3.

To change to install mode for Windows Server 2008 R2 installations:

- 1. Click Start > Run.
- 2. Type change user /install. (Note: There is one space between "user" and "/".)
- 3. Click OK.

Other Steps to Take Before You Begin the Installation

- If you are upgrading from Cornerstone 7.6 and you plan to install on a Windows Vista workstation, you need to uninstall the Cornerstone PDF Printer prior to installing Cornerstone 8.3. To uninstall, go to Control Panel (in Classic View) > Programs and Features > select the Cornerstone PDF Printer and click the Uninstall/Change button.
- If you are using the Mobile Computing module, bring in, synchronize and undo your mobile computers. (See chapter 7, "Installation for the Mobile Computing Module," on page 33 in this guide for more information.)
- 3. Run End of Day processing and perform your backup.
- 4. Close any programs that are open, including Cornerstone, on all computers.
- 5. Close the Cornerstone database engine on the server by right-clicking on the **SQL** icon in the notification area on your taskbar and selecting **Exit**.

Note: To verify the database is closed (in the event the icon is not present), press CTRL+ALT+DELETE and then click **Task Manager**. On the **Processes** tab, select the **dbsrv7.exe** process in the list (for a stand-alone system, select **dbeng7.exe**), click **End Process**, and then click **Yes**.

A message may appear stating that one or more connections are open and verifying that you wish to close the database anyway. If you have verified that Cornerstone has been closed on all computers, click **Yes**.

5: Installation on a Server or Stand-alone Computer

Installing Cornerstone 8.3 from a DVD

The steps in this section cover using a DVD to install a new copy of the Cornerstone* Practice Management System or upgrade Cornerstone in the following installation environments:

- Server
- Remote Desktop Services/Terminal Services server
- Stand-alone computer

If you have more than one computer on your Cornerstone network, install Cornerstone 8.3 on the server first. Then, see chapter 6, "Installation on Workstations," starting on page 28 for information on the proper installation procedures for Cornerstone workstations.

To install Cornerstone from a DVD:

1. Insert the IDEXX Cornerstone DVD in the DVD drive.

If the DVD does not automatically start, click the **Start** button and select **Run**. Type the letter of the DVD drive followed by :\setup.exe (e.g., D:\setup.exe) and click **OK** to open the install program.

Note: If a DVD drive is not available on the server (or stand-alone computer), Cornerstone can be installed across the network from a shared DVD drive on another computer by running the setup.exe file.

- 2. When the Welcome window appears, click Next.
- Click I accept the terms of the license agreement to accept the software license agreement, and then click Next.
- 4. On the Activation Key-Edit Data window, do one of the following:
 - For upgrades: If you are upgrading from an earlier version of Cornerstone, the Clinic ID and Key text boxes will be filled in automatically. Click **Next**.
 - For new installations or installations with a new key: If Cornerstone has never been installed on this computer before, type your clinic ID and key and click **Next**.

Activation K	ey 🛛 🔀
Edit Data Enter requ	ested data.
	clinic ID and activation key, then click Next to continue. A clinic ID and activation obtained by contacting Technical Support at (800) 695-2877.
Clinic ID:	
Key:	l
InstallShield —	
	< Back Next > Cancel

Activation Key—Edit Data window

- 5. If you are installing a multi-computer system, the Setup Type window appears. On the Setup Type window, do one of the following:
 - For upgrades: If you are upgrading from an earlier version of Cornerstone, the installation type will be selected for you automatically. Leave this selection at its auto-detected setting and click **Next**.
 - For new installations: If Cornerstone has never been installed on this computer, select the appropriate installation type and click **Next**.

IDEXX Cornerstone Setup	
Setup Type	124
Please select type of installation.	
O Server	
◯ Workstation	
 Terminal Services Server 	
◯ WinConnect Server	
InstallShield	< Back Next > Cancel

Setup Type window

6. If you have purchased the Central/Remote module, the Setup Type Central window appears. On the Setup Type Central window, select the appropriate Central/Remote installation type and click **Next**.

IDEXX Cornerstone Setup
Setup Type Central
Select Server if you are installing on the file server. Select Workstation if you are installing on an individual workstations. Select Single User if you are installing as a stand alone.
◯ Single User
⊙ Server
O Workstation
InstallShield

Setup Type Central window (appears only for Central/Remote module users)

7. If the IDEXX SmartService* Solutions tool has not been installed with a previous installation of Cornerstone, the IDEXX SmartService License Agreement window appears. (IDEXX SmartService is an Internet-based service and support solution that provides automatic download and notification of new Cornerstone releases and service packs.) Read the license agreement, and then click **Yes** to accept the terms of the agreement and install IDEXX SmartService. A message states that the Cornerstone database will be backed up.

8. Click OK.

A message states that the backup was successful.

9. Click OK.

The Select Features window appears.

10. On the Select Features window, scroll through the list of Cornerstone feature modules and ensure that the appropriate modules are selected (the Cornerstone modules you have purchased are automatically selected). We do not recommend deselecting any module.

IDEXX Cornerstone Setup	
Select Features Select features to install	
Select the features you want to install, and deselect the features you do not want to install.	
Foundation	<u>^</u>
Diagnostic Imaging	
✓ Patient Advisor	≡
☑ In-Clinic Lab Interface	
Practice Explorer	
✓ Time Clock	
🗌 Central	
Destination Folder C:\Cstone Browse	
InstallShield	el

Select Features window for Cornerstone module installation

Note: If you are using the Central/Remote module, you must click the associated check boxes on the Select Features window. The Central/Remote module features are not selected automatically based on the modules you have purchased.

11. Click Next.

The Start Copying Files window appears.

12. Click Next.

Note: If a Windows File Protection message displays, click **Cancel** to avoid restoring the original versions of the files, and then click **Yes**. These files will not affect stability; the files have been tested by IDEXX and are actually newer versions.

Installation begins. A window displays the installation progress.

- 13. If you have purchased the IDEXX Diagnostic Imaging module, a message appears asking if you want to install DICOM[®] services. Do one of the following:
 - If you currently use DICOM services or plan to do so in the future, click Yes.
 - If you do not use DICOM services, click No.
- 14. If you are installing on a server or Remote Desktop Services/Terminal Services server, a message appears to remind you that the Cornerstone database must be updated after you restart the computer. Click **OK**.

15. On the Setup Complete window, click **Yes**, **I want to restart my computer now** and then click **Finish** to restart your computer.

Note: If you have a Remote Desktop Services/Terminal Services server, you will have to restart manually.

16. After restarting your computer, proceed to the "To update your database" section that follows to update your database.

To update your database:

Note: For installations on Windows Vista, Windows 7 or Windows Server 2008 R2, User Account Control (UAC) must be disabled during the update of the database. The update will not start with UAC enabled. See the "Disabling User Account Control" section on page 19 for steps to disable UAC.

- 1. After completing the Cornerstone 8.3 installation and restarting your computer, log on to the server, Remote Desktop Services/Terminal Services server, or stand-alone computer as usual.
- 2. Do the following to start up the database:
 - If you are installing on a Remote Desktop Services/Terminal Services server, start the database manually by double-clicking the **Cornerstone Database** icon **W** on your desktop.
 - If you are installing on a server or standalone computer, the database should start automatically. However, if you do not see the SQL are icon in the notification area on your taskbar, start the database manually by double-clicking the Cornerstone Database Server icon on your desktop.

The IDEXX Database and Services window appears, and displays a list of starting and stopping activities. After a few moments, the Starting Cornerstone Database message appears, along with a second window. This second window closes automatically, and the IDEXX Database and Services window displays several starting messages, followed by the word Finished.

- 3. Click Close (X) in the upper right of the window to close it.
- Double-click the Cornerstone icon on the desktop or select Start > Programs > Cornerstone > Cornerstone to start Cornerstone.

The following message appears: The executable version (8.3.xx) and the database version do not match. Do you want the database brought up to date?

5. Click Yes.

The Cornerstone Update window appears. All of the feature modules you installed are checked.

Note: If your practice uses the IDEXX SmartLink In-house Laboratory module, you must shut down the IDEXX VetLab Station and all other connected analyzers at this time to avoid errors while the database update program is running.

Note: If your images are being stored in a location other than your server, that workstation must be logged into Windows at this time (before updating the database) to avoid altering your Cornerstone default image location.

6. Click **OK** to start the database update process.

Note: Before the database update can begin, Cornerstone must first validate the directories required by Cornerstone. If any validation issues are encountered, a Cornerstone directory validation window appears. Review the information provided and then click **Close** and perform the suggested troubleshooting tasks. Click **Update** if you have completed the troubleshooting tasks and you are

ready to update the Cornerstone database. If you need additional troubleshooting help, please contact IDEXX Cornerstone Support at 1-800-695-2877.

- During the database update, the update program checks to verify the location of multiple folders. If the update detects any changes, you will be prompted multiple times to update the location. Click Yes at each prompt.
- 8. During the database update, the update program checks to verify the type of backup software installed on your server or stand-alone system.
 - If the update program locates one of the supported backup software packages, select the appropriate backup software and click **OK**.
 - If your practice uses other backup software, the software will continue to operate normally (your backup procedure will not change). Click **OK** to continue with the update.
- 9. When the database update is complete, click **OK**.

Note: The first time a staff member logs into Cornerstone after loading Cornerstone 8.3, the Physical Exam Upload window may prompt you to specify the name your practice uses for each of the given species. (If your practice has already established species links, this window does not appear.)

10. In Cornerstone, a Cornerstone administrator will need to adjust security settings for your practice so that staff members can access the new features. See chapter 8, "Security and Setup," for a listing of these security settings.

Upgrading to Cornerstone 8.3 from the Internet using IDEXX SmartService* Solutions

The steps in this section cover using IDEXX SmartService* Solutions to upgrade Cornerstone in the following installation environments:

- Server
- Remote Desktop Services/Terminal Services server
- Stand-alone computer

If you have more than one computer on your Cornerstone network, install Cornerstone 8.3 on the server first. Then, see chapter 6, "Installation on Workstations," starting on page 28 for information on the proper installation procedures for Cornerstone workstations.

The IDEXX SmartService Solutions tool provides automatic downloads and notification of Cornerstone updates. This automated process saves you time, while also ensuring that you are taking advantage of the latest updates and capabilities available in your IDEXX Cornerstone software. The time required to complete the update installation will vary, but in general, the total time required should be less than (or comparable to) the time it typically takes you to install a Cornerstone upgrade.

How it works

When a new Cornerstone update is available, it downloads automatically behind the scenes on the Cornerstone server (or on a stand-alone computer). Once an update is received, all Cornerstone administrators will receive a message upon login to notify them that an update is available and is ready to be installed on the server.

Note: The message appears only once. After closing the message, a yellow **!!** notification icon displays on the Cornerstone status bar (lower right) as a reminder until the update is installed.

When ready to install the update on the server, make sure you are logged in as a Windows®

administrator and double-click the **Update Launcher** icon on the server's desktop to begin the installation.

Once the update has been installed on the server, each Cornerstone workstation will be prompted to update as Cornerstone is opened. The logged-in user can simply click **Yes** to start the Install from Server feature, which will automatically install the update.

To install the Cornerstone update on the server or a stand-alone system:

IMPORTANT: You will need to restart after installation.

Getting Started

- 1. Close any programs that are open, including Cornerstone, on all computers.
- 2. Close the Cornerstone database engine on the server by right-clicking on the **SQL W** icon in the notification area on your taskbar and selecting **Exit**.

Note: To verify the database is closed (in the event the icon is not present), press CTRL+ALT+DELETE and then click **Task Manager**. On the **Processes** tab, select the **dbsrv7.exe** process in the list (for a stand-alone system, select **dbeng7.exe**), click **End Process**, and then click **Yes**.

A message may appear stating that one or more connections are open and verifying that you wish to close the database anyway. Click **Yes**.

Installation

Note: If you have Microsoft[®] Windows Server[®] 2008 R2 installed, you must install the update while in install mode for the installation to work correctly. To change to install mode: 1) Click **Start > Run**. 2) Type **change user /install**. (Note: There is one space between "user" and "/".) 3) Click **OK**.

- 1. Double-click the **Update Launcher I** icon on the desktop.
- 2. On the Cornerstone Update Available dialog box, click Install.
- 3. Follow the prompts to install the Cornerstone update.
- 4. When the installation is complete, click **Finish** to restart the server. This may take a few moments.
- 5. Once the server has been restarted, log in to Windows.

Updating the Database

Before You Begin:

- If your practice uses the IDEXX SmartLink* In-house Laboratory module, you must shut down the IDEXX VetLab* Station and all other connected analyzers at this time to avoid errors while the database update program is running.
- If your images are being stored in a location other than your server, that workstation must be logged in to Windows at this time (before updating the database) to avoid altering your Cornerstone default image location.
- 1. Do the following to start up the database:
 - If you are installing on a Remote Desktop Services/Terminal Services server, start the database manually by double-clicking the Cornerstone Database icon an your desktop.
 - If you are installing on a server or standalone computer, the database should start automatically. However, if you do not see the SQL icon in the notification area on your

taskbar, start the database manually by double-clicking the **Cornerstone Database Server** icon on your desktop.

The IDEXX Database and Services window appears, and displays a list of starting and stopping activities. After a few moments, the Starting Cornerstone Database message appears, along with a second window. This second window closes automatically, and the IDEXX Database and Services window displays several starting messages, followed by the word Finished.

- Click Close (X) in the upper right of the window to close it. Double-click the Cornerstone icon on the desktop or select Start > Programs > Cornerstone > Cornerstone to start Cornerstone.
 The following message appears: "The executable version and the database version do not match. Do you want the database brought up to date? "
- 3. Click **Yes**. The Cornerstone Update window appears. All of the feature modules you installed are checked.
- 4. Click **OK** and follow the prompts through the update process.

6: Installation on Workstations

The steps in this section apply to the following installation environments:

- New Cornerstone workstations (Cornerstone* Practice Management Software has not been installed on this workstation before)
- Upgrades from Cornerstone 7.5 and earlier versions
- Upgrades from Cornerstone 7.6, 8.0, 8.1, or 8.2

The process for upgrading from Cornerstone version 7.6, 8.0, 8.1, or 8.2 begins on this page. The installation process for first-time installations or upgrading from 7.5 and earlier versions begins on page 29.

Upgrades from Cornerstone 7.6, 8.0, 8.1, or 8.2

Starting with Cornerstone version 7.6, when you installed Cornerstone manually on the server and workstations, the Install from Server service was automatically loaded on each machine. That installation put the service in place for future installations so that any future versions/updates/service packs need to be installed only on the server.

Notes on Install from Server:

- The Install from Server service applies to version upgrades only (full release upgrades, minor version updates and service packs). It does **not** apply to new workstation installations of Cornerstone. New installations must be installed from DVD or installed manually (by running the setup.exe file) from a shared network folder.
- The activation key will be acquired automatically from the server.
- On each workstation, the service will use the same install folder location that was used for previous Cornerstone installations on that workstation.

Installing on Workstations Using Install from Server

Note: Before performing Install from Server updates, we recommend that you disable screen savers on all workstations to ensure that the updates can install without interruption.

- 1. Ensure that Cornerstone 8.3 has been installed on the server and that the Cornerstone database has been updated on the server.
- On the workstation, double-click the Cornerstone icon on the desktop or select Start > Programs
 > Cornerstone > Cornerstone to launch Cornerstone.

The following message appears: "The Cornerstone server is running a more recent version of Foundation than this computer. An update was found on the server. Do you wish to install it now?"

3. Click **Yes**. Cornerstone closes and the new version is installed automatically from the server (you do not need to click through the usual installation wizard screens).

Note: If you click **No** when asked if you wish to install the update, Cornerstone closes. Cornerstone will not run on the workstation until the system detects that the workstation version matches the server version.

4. When the installation is complete, restart the workstation.

Installing on Workstations from a Shared Network Folder

If you have problems installing on workstations using the Install from Server service (because of firewall or security issues, etc.), Cornerstone can be installed on the workstations across the network from a shared network folder.

- 1. Ensure that Cornerstone 8.3 has been installed on the server and that the Cornerstone database has been updated on the server.
- 2. On the workstation, browse to the \\server\cstoneshare\install\release shared network folder and double-click the setup.exe icon.

Note: If your server name is something other than *server*, replace the word *server* with your server name. To check the server name: On the server desktop or **Start** menu, right-click the **My Computer** icon and select **Properties**. The server name can be found on the **Computer Name** tab, listed after *Full computer name*.)

- 3. Follow the prompts to install on the workstation.
- 4. When the installation is complete, restart the workstation.

New Workstations or Upgrades from Cornerstone 7.5 and Earlier Versions

1. Insert the IDEXX Cornerstone DVD in the DVD drive.

If the DVD does not automatically start, click the **Start** button and select **Run**. Type the letter of the DVD drive followed by :\setup.exe (e.g., d:\setup.exe) and click **OK** to open the install program.

Note: If a DVD drive is not available on this workstation, Cornerstone can be installed across the network from a shared DVD drive on another computer by running the setup.exe file.

- 2. When the Welcome window appears, click Next.
- 3. Click **I accept the terms of the license agreement** to accept the software license agreement, and then click **Next**.
- 4. On the Activation Key-Edit Data window, do one of the following:
 - For upgrades: If you are upgrading from an earlier version of Cornerstone, the Clinic ID and Key text boxes will be filled in automatically. Click **Next**.
 - For new installations or installations with a new key: If Cornerstone has never been installed on this computer before, type your clinic ID and key and click **Next**.

Note: If you do not know the clinic ID or key, you can acquire it from the server. To access this information, log in to Cornerstone on the server, then go to **Help > About IDEXX Cornerstone** and click the **More** button.

Activation K	ey 🛛 🛛
Edit Data Enterrequi	eter haise
	clinic ID and activation key, then click Next to continue. A clinic ID and activation obtained by contacting Technical Support at (800) 695-2877.
Clinic ID:	
Key:	
InstallShield —	< Back Next > Cancel

Activation Key—Edit Data window

5. On the Setup Type window, if not already selected, select the **Workstation** installation type and click **Next**.

IDEXX Cornerstone Setup	$\overline{\mathbf{X}}$
Setup Type	
Please select type of installation.	
O Server	
 Workstation 	
 Terminal Services Server 	
◯ WinConnect Server	
InstallShield	
	< Back Next > Cancel

Setup Type window

6. If you have purchased the Central/Remote module, the Setup Type Central window appears. On the Setup Type Central window, select the **Workstation** installation type and click **Next**.

IDEXX Cornerstone Setup	
Setup Type Central	
Select Server if you are installing on the file server, an individual workstations. Select Single User if yo	
🔿 Single User	
◯ Server	
 Workstation 	
InstallShield	<back next=""> Cancel</back>

Setup Type Central window (appears only for Central/Remote module users)

The Select Features window appears.

7. On the Select Features window, scroll through the list of Cornerstone feature modules and ensure that the appropriate modules are selected (the Cornerstone modules you have purchased are automatically selected). We do not recommend deselecting any module.

IDEXX Cornerstone Setup	×
Select Features Select features to install	
Select the features you want to install, and dea	select the features you do not want to install.
Foundation	<u>^</u>
Diagnostic Imaging	
Patient Advisor	
In-Clinic Lab Interface	
Practice Explorer	
Time Clock	
Central	
	~
Destination Folder	
C:\Cstone	Browse
InstallShield	
	< Back Next > Cancel

Select Features window for Cornerstone module installation

Note: If you are using the Central/Remote module, you must click the associated check boxes on the Select Features window. The Central/Remote module features are not selected automatically based on the modules you have purchased.

- 8. Click Next.
- 9. The Start Copying Files window appears. Click Next.

Note: If a Windows File Protection message is displayed, click **Cancel** to the restoration of the original versions of these files and then click **Yes** that you are sure. These files will not affect stability; the files have been tested by IDEXX and are actually newer versions.

Installation begins. A window displays the installation progress.

10. On the Setup Complete window, click **Yes**, **I want to restart my computer now** and then click **Finish** to restart your computer.

7: Installation for the Mobile Computing Module

Installing Cornerstone 8.3 for the Mobile Computing Module

If you are using the Mobile Computing module, use the following instructions to install the Cornerstone* 8.3 Practice Management System.

Note: See the *Cornerstone User's Manual* for additional information on using the Mobile Computing module.

- 1. Bring all mobile computers into the practice. (If you cannot bring all mobile computers into the practice at the same time, see the "If You Cannot Bring Each Laptop/Tablet into the Practice" section that follows.)
- 2. Synchronize each mobile computer.
- 3. Remove each mobile computer from the Mobile Computing List by selecting **Undo mobile computer**.
- 4. At your practice, install Cornerstone 8.3 on your server and workstations according to the steps in chapter 5, "Installation on a Server or a Stand-alone Computer," and chapter 6, "Installation on Workstations," in this guide.
- 5. Install Cornerstone 8.3 on each laptop/tablet according to the steps in chapter 6, "Installation on Workstations," starting on page 28 in this guide.
- 6. Start Cornerstone on each laptop/tablet and create the mobile computer. (See the "Using the Mobile Computing Module" chapter in the *Cornerstone User's Manual* for information on creating a mobile computer in Cornerstone.)

If You Cannot Bring Each Laptop/Tablet into the Practice

We strongly encourage you to bring all of the mobile computers into the practice at the same time to perform the Cornerstone 8.3 installation.

WARNING: If you do not bring in all mobile computers to be synchronized at the time of updating, you will be at risk of data loss on the mobile computers that were not present at the time of the update.

If it is not possible to bring all mobile computers into the practice at the same time, follow these instructions:

- 1. Bring in as many of your mobile computers as possible and follow the instructions in the previous "Installing Cornerstone 8.3 for the Mobile Computing Module" section.
- 2. As you are able, bring your other mobile computers into the practice for updating.

You will need to update your mobile database before synchronizing:

- a. On each mobile computer that has not yet been updated, install Cornerstone 8.3 according to the steps in chapter 6, "Installation on Workstations," starting on page 28 in this guide.
- b. After installation, start Cornerstone on the mobile computer. An alert will appear, informing you that your database needs to be updated. Click **Yes** to update the database. (If the update is not performed, Cornerstone will not function on this mobile computer and you will not be able to synchronize with the server.)

- 3. Synchronize the mobile computer with the server.
- 4. For any mobile computer that is not synchronized with the server before Cornerstone 8.3 is installed on the server, the first time you sync after installing 8.3 on the mobile computer, you will see a duplicate client error for EMR_CLT and a duplicate patient error for EMR_PAT. To solve this problem, do the following:
 - a. In order to complete the synchronization, enter a new patient ID and a new client ID (you will need to delete these IDs in the next step, so you may want to write them down for easy reference).
 - Log in to Cornerstone on the server and delete the new patient ID and client ID (using Lists > Patient and Lists > Client). Note that the patient ID must be deleted first, as you cannot delete a client ID that has a patient ID associated with it.

Note: Do not delete the EMR CLT or EMR PAT IDs from the mobile computer or the server.

8: Security and Setup

Suggested Security by Staff and Role for New Windows in Cornerstone 8.3

The table below shows the suggested security access for different staff and roles within your practice:

- **GEN**—If you choose to have only two tiers of security, for administrators and all other staff members, this would be the access settings for staff members who aren't administrators.
- **REC**—Receptionists
- **DR**—Doctors
- **TECH**—Technicians
- M/O— Managers/Owners

For more information about setting up security, see the Cornerstone 8.3 Default and Security Settings Guide.

Dialog Security		Acce	ss Sug	gestions		Path to Locate
	GEN	REC	DR	TECH	M/O	
Antech – Lab Item Price Update			Х	Х	Х	Activities > Lab Work > Update Antech Lab Item Prices
Imaging DICOM Modality Types			Х	Х	Х	Activities > Imaging > DICOM Modality Types
Invoice Discount – Apply to all lines	х	х	х	х	х	From the Invoice, select No in Discount column > Select discount to be applied > Select option to Add line discounts to all lines for this patient > Select OK
IDEXX Message Center			Х	Х	Х	፤≡ button on toolbar
Patient Vaccine Corrections			Х	Х	х	Select patient Tags/Microchip ID tab > right-click on vaccination > Correct or Reissue
Synchronize LifeLearn Documents					Х	Tools > Synchronize LifeLearn Documents
Report Security		Acce	ss Sug	gestions		Path to Locate
	GEN	REC	DR	TECH	M/O	
EOM Referral Revenue Report by Doctor by Revenue Center					Х	Reports > End of Period > EOM > Specify Monthly Period>Specify EOM - Referral Revenue Report by Doctor by Revenue Center> Print OR Preview>Specify Doctor/Doctors>OK
EOM Referral Revenue Report by Hospital by Revenue Center					х	Reports > End of Period > EOM > Specify Monthly Period > Specify EOM - Referral Revenue Report by Hospital by Revenue Center > Print OR Preview > Specify Hospital or Revenue Center
EOM Referral Revenue Report by Rev Center by Doctor by Hospital					х	Reports > End of Period > EOM > Specify Monthly Period > Specify EOM Referral Revenue Report by Rev Center

			by Doctor by Hospital > Print OR Preview > Specify Date, Doctor ID, or Revenue Center ID
EOM Referral Revenue Report by Rev Center by Hospital by Doctor		х	Reports > End of Period > EOM > Specify Monthly Period>Specify EOM Referral Revenue Report by Rev Center by Hospital by Doctor>Specify Hospital ID or Revenue Center ID
EOM Referral Summary by Revenue Center – Doctor by Hospital		Х	Reports > End of Period > EOM > Specify Monthly Period>Specify EOM Referral Summary by Revenue Center – Doctor by Hospital> Print or Preview>Specify Doctor ID or Revenue Center ID
EOM Referral Summary by Revenue Center – Hospital by Doctor		Х	Reports > End of Period > EOM > Specify Monthly Period>Specify EOM Referral Summary by Revenue Center – Hospital by Doctor>Print or Preview>Specify Hospital ID or Revenue Center ID

This document contains suggested security settings; however, it is ultimately the practice's decision-maker that decides who needs access to certain areas within Cornerstone, and IDEXX shall not be liable for any actions arising from those security decisions.

9: Communications and Training

Training Checklist

Use the "Use" sections in the Cornerstone* 8.3 What's New to train your staff on the new features.

The table below provides a checklist for all the new features in What's New.

Image: Properties of the second client of the second client response of t	Client Communications					
IDEXX Pet Health Network* Pro Automatic Appointment Confirmations through IDEXX Pet Health Network Pro Share Documents with Clients New Features to Record Client Email Addresses IVewing IDEXX Pet Health Network Pro Activities on Reports Communications Management for Email and Fax Overview Reminder Emailing or Faxing Documents from the Cornerstone Editor Viewing Communications Log Accessing the Communications Log Next Appointment Scheduling Reminder Scheduling an appointment using the reminder Additional Client Information Enhancements New Features to Help Record Client Email Addresses Infancements Automatically Printing a Linked "Diagnosis Document" from a Medical Note Diagnosis Codes Customizing and Editing Your Library of ClientEd Online Documents Infancements Selecting a LifeLearn Document from the Start New Document Window IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test IDEXX VetLab* Station Verview IDEXX VetLab Verview IDEXX VetLab Verview IDEXX Pet Health Network Pro Verview IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test						
Pro Share Documents with Clients New Features to Record Client Email Addresses Viewing IDEXX Pet Health Network Pro Activities on Reports Communications Management for Email and Fax Overview Emailing or Faxing Documents from the Cornerstone Editor Viewing Communication in the Patient Clipboard Viewing the Communications Log Accessing the Communications Shows staff who entered confirmation Appointment confirmations shows staff who entered confirmation Automatically Printing a Linked "Diagnosis Document" from a Medical Note Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Online Documents Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Diter Ways to Use LifeLearn ClientEd Documents with Cornerstone Selec		View Incoming Client Requests and information via the IDEXX Message Center				
Communications Management Image Decembers on Record Client Email Addresses Communications Management Image Decembers on Record Client Email Addresses Communications Management Image Decembers on the Connerstone Editor for Email and Fax Image Decembers on the Connerstone Editor Viewing Communication in the Patient Clipboard Viewing Communications Log Next Appointment Scheduling Image Decembers on Scheduling an appointment using the reminder Additional Client Information Image Decembers on Scheduling an appointment using the reminder Additional Client Information Image Decembers on Scheduling an appointment using the reminder Linking Documents to Diagnosis Codes Image Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online With Connerstone Image Decembers on Scheduling Your Library of ClientEd Online Documents Integration with Cornerstone Image Decembers on Scheduling Your Library of ClientEd Online Documents Integration with Cornerstone Image Decembers on Scheduling Your Library of ClientEd Online Window Integration with Cornerstone Image Decembers on Scheduling Your Library of ClientEd Online Documents Integration with Cornerstone Image Decembers on Scheduling Your Library of ClientEd Online Documents Integration with Cornerstone Image Decembers on D	IDEXX Pet Health Network*	Automatic Appointment Confirmations through IDEXX Pet Health Network Pro				
Image: constraint of the second sec	Pro	□ Share Documents with Clients				
Communications Management for Email and Fax Overview Emailing or Faxing Documents from the Cornerstone Editor Viewing Communication in the Patient Clipboard Viewing the Communications Log Accessing the Communications Log Next Appointment Scheduling Reminder Cheduling an appointment using the reminder Additional Client Information New Features to Help Record Client Email Addresses Appointment confirmations shows staff who entered confirmation Linking Documents to Diagnosis Codes Automatically Printing a Linked "Diagnosis Document" from a Medical Note Customizing and Editing Your Library of ClientEd Online Documents Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory VetConnect* PLUS Coverview Understanding the VetConnect PLUS Tab Dielect an IDEXX VetLab Station when requesting a test Overview Viewing and Updating Profiles Ordering a Profile Updating a Pending Test Request 		New Features to Record Client Email Addresses				
Communications Management for Email and Fax Emailing or Faxing Documents from the Cornerstone Editor Viewing Communication in the Patient Clipboard Viewing the Communications Log Accessing the Communications Log Accessing the Communications Log Rest Appointment Scheduling meminder Scheduling an appointment using the reminder Additional Client Information Enhancements Appointment confirmations shows staff who entered confirmation Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online^w Integration with Cornerstone Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Understanding the VetConnect PLUS Tab Inderstanding the VetConnect PLUS Tab IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test Overview Viewing and Updating Profiles Ordering a Profile Updating a Pending Test Request		Viewing IDEXX Pet Health Network Pro Activities on Reports				
Communications Management for Email and FaxUsewing Communication in the Patient Clipboard Usewing the Communications Log Accessing the Communications LogNext Appointment Scheduling ReminderScheduling an appointment using the reminderAdditional Client Information EnhancementsNew Features to Help Record Client Email Addresses Appointment confirmations shows staff who entered confirmationLinking Documents to Diagnosis CodesAutomatically Printing a Linked "Diagnosis Document" from a Medical NoteLifeLearn® ClientEd Online™ Integration with CornerstoneCustomizing and Editing Your Library of ClientEd Online Documents Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with CornerstoneVetConnect* PLUSOverview Understanding the VetConnect PLUS TabIDEXX VetLab* StationSelect an IDEXX VetLab Station when requesting a test Uverview Uverview Overview 						
for Email and Fax I Viewing Communication in the Patient Clipboard I Viewing the Communications Log I Accessing the Communications Log Next Appointment Schedulling I Scheduling an appointment using the reminder Additional Client Information I New Features to Help Record Client Email Addresses Enhancements I Automatically Printing a Linked "Diagnosis Document" from a Medical Note Linking Documents to I Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online™ I Customizing and Editing Your Library of ClientEd Online Documents Integration with Cornerstone I Selecting a LifeLearn Document from the Start New Document Window Integration with Cornerstone I Overview VetConnect* PLUS I Overview IDEXX VetLab* Station I Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Viewing and Updating Profiles I Overview I Overview	Communications Management					
Image: Constraint of the communications LogNext Appointment Scheduling ReminderScheduling an appointment using the reminderAdditional Client Information EnhancementsNew Features to Help Record Client Email Addresses Appointment confirmations shows staff who entered confirmationLinking Documents to Diagnosis CodesAutomatically Printing a Linked "Diagnosis Document" from a Medical NoteLifeLearn® ClientEd Online™ Integration with CornerstoneCustomizing and Editing Your Library of ClientEd Online Documents Selecting a LifeLearn Document from the Start New Document Window Doter Ways to Use LifeLearn ClientEd Documents with CornerstoneLaboratoryOverview Understanding the VetConnect PLUS TabVetConnect* PLUSSelect an IDEXX VetLab Station when requesting a testIDEXX VetLab* StationSelect an IDEXX VetLab Station when requesting a testAntech Diagnostics® Module Updating a Profile Updating a Profile Updating a Pending Test Request	_					
Next Appointment Scheduling Reminder Scheduling an appointment using the reminder Additional Client Information Enhancements New Features to Help Record Client Email Addresses Linking Documents to Diagnosis Codes Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online™ Integration with Cornerstone Customizing and Editing Your Library of ClientEd Online Documents LifeLearn® ClientEd Online™ Integration with Cornerstone Customizing and Editing Your Library of ClientEd Document Window Laboratory Other Ways to Use LifeLearn ClientEd Documents with Cornerstone LiDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Overview Updating a Profile Ordering a Profile Updating a Pending Test Request Updating Test Request						
Reminder I Scheduling an appointment using the reminder Additional Client Information Enhancements I New Features to Help Record Client Email Addresses I Appointment confirmations shows staff who entered confirmation Appointment confirmations shows staff who entered confirmation Linking Documents to Diagnosis Codes I Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online" I Customizing and Editing Your Library of ClientEd Online Documents I Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory I Overview Understanding the VetConnect PLUS Tab IDEXX VetLab* Station I Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Overview I Overview I Overview I Overview I Overview I Dietring and Updating Profiles Ordering a Profile I Ordering a Profile I Updating a Pending Test Request		Accessing the Communications Log				
Enhancements□ Appointment confirmations shows staff who entered confirmationLinking Documents to Diagnosis Codes□ Automatically Printing a Linked "Diagnosis Document" from a Medical NoteLifeLearn® ClientEd Online™ Integration with Cornerstone□ Customizing and Editing Your Library of ClientEd Online Documents □ Selecting a LifeLearn Document from the Start New Document Window □ Other Ways to Use LifeLearn ClientEd Documents with CornerstoneLaboratory□ Overview □ Understanding the VetConnect PLUS TabVetConnect* PLUS□ Select an IDEXX VetLab Station when requesting a testIDEXX VetLab* Station□ Overview □ Viewing and Updating Profiles □ Ordering a Profile □ Updating a Pending Test Request		□ Scheduling an appointment using the reminder				
Linking Documents to Image: Automatically Printing a Linked "Diagnosis Document" from a Medical Note LifeLearn® ClientEd Online™ Image: ClientEd Online™ Integration with Cornerstone Image: Selecting a LifeLearn Document from the Start New Document Window Image: Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Image: Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory Image: Overview Image: Overview VetConnect* PLUS Image: Overview Image: Overview IDEXX VetLab* Station Image: Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Image: Overview Image: Overview Image: Overview Image: Overview	Additional Client Information	New Features to Help Record Client Email Addresses				
Diagnosis Codes Automatically Printing a Linked "Diagnosis Document" from a Medical Note Diagnosis Codes LifeLearn® ClientEd Online T Selecting a LifeLearn Document from the Start New Document Window Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory Overview Understanding the VetConnect PLUS Tab IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test Overview Viewing and Updating Profiles Ordering a Profile Updating a Pending Test Request 	Enhancements	Appointment confirmations shows staff who entered confirmation				
LifeLearn® ClientEd Online™ □ Selecting a LifeLearn Document from the Start New Document Window □ Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory VetConnect* PLUS □ Overview □ Understanding the VetConnect PLUS Tab IDEXX VetLab* Station □ Select an IDEXX VetLab Station when requesting a test □ Overview □ Updating a Pending Test Request	0	□ Automatically Printing a Linked "Diagnosis Document" from a Medical Note				
Integration with Cornerstone Selecting a LifeLearn Document from the Start New Document Window Integration with Cornerstone Other Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory Integration of the VetConnect PLUS Tab VetConnect* PLUS Inderstanding the VetConnect PLUS Tab IDEXX VetLab* Station Select an IDEXX VetLab Station when requesting a test Integration Diagnostics® Module Viewing and Updating Profiles Integration a Profile Updating a Pending Test Request		Customizing and Editing Your Library of ClientEd Online Documents				
Image: Constant of the Ways to Use LifeLearn ClientEd Documents with Cornerstone Laboratory VetConnect* PLUS Image: Document of the VetConnect PLUS Tab IDEXX VetLab* Station Image: Document of the VetConnect PLUS Tab IDEXX VetLab* Station Image: Document of the VetConnect PLUS Tab Image: Document of the VetConnect PLUS Tab IDEXX VetLab* Station Image: Document of the VetConnect PLUS Tab Image: Document of the VetConnect PLUS Tab Image: Document of the VetConnect PLUS Tab Image: Document of the VetConnect of the VetConnect PLUS Tab Image: Document of the VetConnect of the VetConn		Selecting a LifeLearn Document from the Start New Document Window				
VetConnect* PLUS Image: Overview Image: Understanding the VetConnect PLUS Tab IDEXX VetLab* Station Image: Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Image: Overview Image: Overview Image: Overview Image: Ove		Other Ways to Use LifeLearn ClientEd Documents with Cornerstone				
VetConnect* PLUS Image: Understanding the VetConnect PLUS Tab IDEXX VetLab* Station Image: Select an IDEXX VetLab Station when requesting a test Antech Diagnostics® Module Image: Overview Image:	Laboratory					
IDEXX VetLab* Station Import Select an IDEXX VetLab Station when requesting a test IDEXX VetLab* Station Import Overview Antech Diagnostics® Module Import Overview Import Overview Import Overview Imp	VetConnect* PLUS					
Antech Diagnostics® Module Overview	Veloonneci FLOS	Understanding the VetConnect PLUS Tab				
Antech Diagnostics [®] Module Uviewing and Updating Profiles Crdering a Profile Updating a Pending Test Request	IDEXX VetLab* Station	□ Select an IDEXX VetLab Station when requesting a test				
Antech Diagnostics® Module □ Ordering a Profile □ Updating a Pending Test Request						
 Ordering a Profile Updating a Pending Test Request 	Antech Diagnostics® Module	Viewing and Updating Profiles				
		□ Ordering a Profile				
Specifying Staff IDs for Laboratory Requests		Updating a Pending Test Request				
		□ Specifying Staff IDs for Laboratory Requests				
Miscellaneous Laboratory Development Result Download and Result Matching	Miscellaneous Laboratory	Result Download and Result Matching				
Enhancements DIDEXX Lab Report	Enhancements	□ IDEXX Lab Report				
		Linking Cornerstone Values with Laboratory Values				
		Linking Cornerstone Values with Laboratory Values				

Imaging		
Cornerstone Imaging	Enhanced Crop Tool	
	Imaging Connections Status Bar	
IDEXX Diagnostic Imaging	□ Improved DICOM [®] Queue	
Diagnostic and IDEXX Digital Radiography	Changes to the Image Capture Window (Overview)	
	□ Image Coach* Tools	
	More Powerful IDEXX I-Vision Mobile* Application	
Modality Worklist Integration		
	Option to Select a Modality When Creating a New Image Request	
	Enhancements to the Image Requests and Results Window	
	Resolving Orphan or Not Requested Statuses for a Non-IDEXX Imaging Modality	
	Manually Setting Requests to a Status of Completed	
	Image Request Special Action Item Report Filtered by Modality	
Credit Card Processing		
Electronic Signature Capture with Integrated Payment Solution	Capturing the Client Signature at Payment Time	
	Posting an Invoice—New Option to Print a Full Invoice or Invoice Receipt	
	□ Identifying Signed Credit Card Transactions on the Patient Clipboard and the Client	
	Account	
	Changes to Payment on Account Report	
Client Communications		
Communications Management for Email and Fax		
	Emailing or Faxing Documents from the Cornerstone Editor	
	□ Viewing Communication in the Patient Clipboard	
	□ Viewing the Communications Log	
	Accessing the Communications Log	
Linking Documents to Diagnosis Codes	Automatically Printing a Linked "Diagnosis Document" from a Medical Note	
LifeLearn [®] ClientEd Online™ Integration with Cornerstone	Customizing and Editing Your Library of ClientEd Online Documents	
	Selecting a LifeLearn Document from the Start New Document Window	
	□ Other Ways to Use LifeLearn ClientEd Documents with Cornerstone	
Referral Relationship Management		
Referral Relationship Management		
	□ Finding and Viewing Referral Veterinarian Information on the Doctor List	
	Adding, Editing and Viewing Referral Veterinarian Information on the Doctor Information Window	
	Finding and Viewing Referral Hospital Information on the Hospital List	
	Adding, Editing, and Viewing Referral Hospital Information on the Hospital Information Window	
	Accessing Referral Information and Referral History in Patient Information	
	Accessing Referral Information on the Patient Clipboard* Window	
	Referral Information on the Patient Check-in/out Window	
	Referral Information for Mobile Computing	
	New and Enhanced Referral Reports	

Inventory	
Enhancements to Cornerstone Inventory Windows	□ Want List
	Purchasing Work List
	Purchase Order List
	Create Purchase Order
	□ Item Information Tabs
	□ Search for Inventory Item
	New Invoice Item Grid Sort and Display Options
Electronic Inventory Ordering with IDEXX SmartOrder* Inventory Solution and MWI® Veterinary Supply	Creating an MWI Purchase Order in Cornerstone
	Adding Vendor Item IDs for Matching Inventory Items in Cornerstone
	Adding a Vendor Item ID While Working with Inventory Orders
	Tracking Order Statuses from Cornerstone
	Synchronization of Electronic Inventory Orders with MWI
Miscellaneous	Login/Password Management
	New Database and Services Management Window
	Additional Changes and Enhancements

Training Files

Work with the Cornerstone tutorial training files data to complete exercises that mimic everyday tasks in your practice without jeopardizing the integrity of your real medical records. You can access the training files from the Microsoft[®] Windows start menu. Select **Start > Programs > Cornerstone > Training Files**. Cornerstone opens, and uses the training files database rather than the live version of the database that contains your practice's information.

In the training files (including the Cornerstone Assistant training files), functionality is limited for some Cornerstone 8.3 features which ordinarily communicate with systems outside of Cornerstone:

- Pet Health Network Pro—You can see and use Pet Health Network Pro features like the share button and the IDEXX Message Center alerts for appointment requests; however, these actions won't be reflected on pethealthnetworkpro.com or petly.com.
- Inventory Ordering using IDEXX SmartOrder with MWI Veterinary Supply—Although you can't use IDEXX SmartOrder, you can view the windows used for working with this feature. Also, Purchase Order 1000 in the training files shows an example of a purchase order created using IDEXX SmartOrder.
- **IDEXX Faxing Service**—Although you can't send faxes from the training files with the IDEXX Faxing Service, this feature is very close to the new email feature that is usable in the training files. Also, completed faxes are shown on the Communications Log in the training files.
- Antech Diagnostics Module—A training version of this module is not available in the training files. However, training is included when you purchase the Antech module or upgrade from the previous version of the module.
- Lifelearn ClientEd Online—You can only set up and synchronize Cornerstone with LifeLearn ClientEd Online documents from your live Cornerstone database.

- VetConnect PLUS—You can view VetConnect PLUS information for the following patients:
 - o Zeke—Patient ID 15329 (Feline)
 - o Blake—Patient ID 3679K-1 (Canine)
 - o Twister-Patient ID 120Q-1 (Canine)
- Electronic signature capture for the integrated payment solution with XCharge[®]—In the training files, you can view XCharge features except the ability to capture electronic signatures.
- Activation—You cannot activate or set up any of the above features from the training files, only from your live Cornerstone database. Also, in the configuration windows for a multi-location single database environment, the training files do not allow you to apply primary settings to a secondary practice.

Appendix: Starting the Database

Starting the Database

The database can only be opened or closed from the server. The database is the only thing required to be running on the server; Cornerstone does not need to be running. If the database is not open, workstations will not be able to access Cornerstone.

The Cornerstone database is the storage location for your data and is housed on the Cornerstone server.

Initially the database is populated with information that is provided during the conversion process from your previous software when you convert to Cornerstone. If you were not previously computerized, more manual entry of data may be required.

The database will continue to be populated by information that you enter during your daily activities.

Generally, Cornerstone's database will automatically start (open) when the server is started. Depending on the type of system setup you have, there are different ways to manually start the Cornerstone database on your server.

Depending on the type of system setup you have, there are different ways to manually start the Cornerstone database on your server.

You will know the database has started when you see the Database Icon (SQL) on the task bar next to the time display in the lower right corner of your screen.

You will know the database has started when you see the SQL icon in the notification area on your taskbar:



You will see this database icon if you have a Peer-to-Peer, Dedicated or Terminal Server.

You will see this database icon if you have a Stand-alone Computer.

Automatic Start

Cornerstone's database will automatically start (open) when the server is started.

The IDEXX Database and Services window appears, and displays a list of starting and stopping activities. After a few moments, the Starting Cornerstone Database message appears, along with a second window. This second window closes automatically, and the IDEXX Database and Services window displays several starting messages, followed by the word Finished.

If the window does not close automatically, click **Close (X)** in the upper right of the window to close it. Only do this after Finished is displayed.

Manual Start

- 1. Do the following to start up the database:
 - On a Remote Desktop Services/Terminal Services server, start the database manually by double-clicking the **Cornerstone Database** icon **Services** on your desktop.

- On other servers or standalone computers, the database should start automatically. However, if you do not see the **SQL** are icon in the notification area on your taskbar, start the database manually by double-clicking the **Cornerstone Database Server** icon on your desktop.
- In either case, if the icon is not available, click **Start** in the lower left of the screen, then select **Programs**>**Cornerstone**> **Cornerstone Database Server**

The IDEXX Database and Services window appears, and displays a list of starting and stopping activities. After a few moments, the Starting Cornerstone Database message appears, along with a second window. This second window closes automatically, and the IDEXX Database and Services window displays several starting messages, followed by the word Finished.

2. Click **Close (X)** in the upper right of the window to close it. Only do this after Finished is displayed.