

Cornerstone 8.5 NEXT Software and IDEXX-PACS

4.2 Imaging Software upgrade checklist

Practice name: _____ Number of seats _____
Current Cornerstone® Software version: _____

Plan

- Select your upgrade champion. Name: _____
- Schedule your upgrade. Date: _____ Time: _____
Server: Plan on 1 to 3 hours to upgrade (varies based on database size, current software version, and diagnostic imaging history).
Workstations and capture station: Plan on 30 minutes to 1 hour each. Upgrade can occur on three workstations at once. You must upgrade Cornerstone Software first on each workstation, and then upgrade IDEXX-PACS™ on each workstation used to view diagnostic images.
- Watch the Cornerstone 8.5 NEXT with IDEXX-PACS 4.2 installation tutorial
- Read the [Cornerstone 8.5 NEXT Software Installation Guide](#)
- Read the release notes for the version of Cornerstone you're on now
 - [Cornerstone 8.2 and 8.3 Release Notes](#)
 - [Cornerstone 8.3 NEXT Release Notes](#)
- Confirm hardware and internet browser meet [system requirements](#)
- Educate your staff on the newest features and enhancements in Cornerstone and IDEXX-PACS Software
 - [Cornerstone 8.5 NEXT podcast](#)
 - Snippets: Cornerstone 8.5 NEXT and Digital Imaging enhancements
 - Snippets: Cornerstone 8.5 NEXT and IDEXX Reference Laboratories enhancements

Prepare (day of upgrade)

- Sync and Undo mobile computers (if applicable)
- Complete end-of-day processing (only applicable if upgrading at the end of the day)
- Perform Cornerstone Software backup and image backup
- Close Cornerstone and restart all computers including your server

.....
To access the installation guide, release notes, system requirements, and staff training, visit idexx.com/upgrade3
.....

Upgrade

- Upgrade the Cornerstone server first
- Upgrade Cornerstone workstations
- For workstations where you want image-viewing capability, once you've upgraded Cornerstone, double-click the IDEXX-PACS icon on the desktop to initiate the IDEXX-PACS upgrade

If you need additional help, contact IDEXX Cornerstone Customer Support at **1-800-695-2877**.

For help with IDEXX digital imaging, contact IDEXX Digital Imaging Customer Support at **1-877-433-9948**.